

BLOOD GLUCOSE TESTING

Candidate: _____ Examiner: _____

Pass

Date: _____ Signature: _____

Fail

	Points Possible	Points Awarded
Time Start: _____ Time Limit: 5 min		
Takes, or verbalizes, body substance isolation precautions	1	
Select appropriate materials:	Glucometer	
	Test Strips	
	Lancet	
	Lancet Holder	
	Tissue	
Wash patient's hands with warm water and dry thoroughly. If you use alcohol, let dry.	1	
Place lancet into lancet device	1	
Cock lancet device	1	
Turn glucometer on	1	
Insert test strip	1	
Choose a spot on side of finger	1	
Apply lancet	1	
Apply drop of blood to test strip	1	
Apply tissue to wound site	1	
Dispose of lancet into sharps container	1	
Read and record results	1	
Dispose of strip in appropriate container	1	
Time End: _____ Minimum Score: 15 Total:	18	

Critical Criteria

- _____ Did not take, or verbalize, body substance isolation precautions
- _____ Failure to wash hands or cleanse site
- _____ Failure to obtain enough blood for testing
- _____ Failure to place lancet in sharps container
- _____ Failure to place test strip in appropriate container